<b>Produktname:</b> Silosept®	 Stanze (druckt nicht)	<input type="checkbox"/>	<b>Kunde:</b> schülke	 RAKO ETIKETTEN GMBH & CO. KG MÖLLNER LANDSTRASSE 15 22969 WITZHAVE Tel. +49 (0) 4104 693 0 Fax +49 (0) 4104 1721 E-Mail: info@rako-etiketten.com
<b>Mat.-nr.:</b> 20000097	 schwarz	<input type="checkbox"/>	<b>Etik.-nr.:</b> 305049128_00	
<b>Index:</b> A	 Pantone 300	<input type="checkbox"/>	<b>Version (intern):</b> 03	
<b>Format:</b> 280 x 105 mm	 Pantone 353	<input type="checkbox"/>	<b>Datum:</b> 11.02.15	
<b>Gebinde:</b> 800 g / HE-DS	 Pantone 032	<input type="checkbox"/>	<b>Grafiker:</b> JR (.indd CS6)	
<b>Land:</b> DE, FR, EN			<b>Schrift &lt; 6pt:</b> 5,5 pt	

# Silosept®

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Immersion disinfectant for dental impressions  
and restorations.

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For professional use only. Caution:  
federal (U.S.A.) Law restricts this  
device to sale by or on the order  
of a dentist.

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## Immersion disinfectant for dental impressions and restorations.

**Range of indications:** Immersing silicone, vinylsiloxanether, polyether and alginate dental impressions as well as dental restorations to disinfect. **Properties:** The mixed solution is microbiologically highly effective and cleans well in only 10 minutes immersion time – it is compatible with impression materials, which remain dimensionally stable. **Mode of application:** Sprinkle 2 scoop (20 g) of Silosept® into 2 litres of lukewarm water (creates a 2 % solution). Rinse the impressions and restorations in running water and immerse them in the solution. Allow 10 minutes for it to take effect, remove them and rinse again in running water. **Change disinfection solution daily. Microbiological effectiveness:** A 2 % Silosept® solution deactivates the relevant germ spectrum – bacteria (incl. Tb) and fungi, enveloped viruses\* (incl. HIV, HBV, HCV) – in 10 minutes, antiviral agent at 2 % - 30 minutes and is easily biodegradable. **Special instructions:** If the patient has been diagnosed with TB, extend the reaction time to 30 minutes and use impression materials that are appropriate for this application. Test material compatibility prior to use. Direct contact with the skin and fabrics may result in bleaching effects. The possibility of sensitive individuals experiencing allergic reactions to the product cannot be ruled out. Wear protective gloves when placing/removing the impressions and restorations. **Composition:** 100 g Silosept® contain: 45 g pentapotassium-bis-(peroxymonosulfate)bis(sulphate) **Labelling according to Regulation (EC) No. 648/2004:** 5 – 15 % anionic surfactants, < 5 % non-ionic surfactants, < 5 % soap and phosphonates, perfumes. **Hazard notes and safety advice: DANGER: H314** Causes severe skin burns and eye damage. **P280** Wear protective gloves (e.g. butyl rubber) / eye protection/face protection. **P301+P310+P330** IF SWALLOWED: Immediately call a POISON CENTER or doctor/ physician. Rinse mouth. **P303+P361+P353** IF ON SKIN (or hair): Remove/ Take off immediately all contaminated clothing. Rinse skin with water/ shower. **P305+P351+P338+P310** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor/ physician. **P501** Dispose of contents/ container to an approved waste disposal plant. Do not store above room temperature and keep closed. Do not use after expiration date.  
\* In accordance with the Robert- Koch- Institute guideline - German Federal Health Gazette 08/2008  
Updated: 10th February 2015

CE 0123  800 g REF 18701

Made in Germany

20000097-A  
LOT  

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**Kettenbach GmbH & Co. KG**