
Sugi®

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

Composition

Cellulose approx.: 35 % (m/m)
 Regenerated cellulose approx.: 65 % (m/m)

The composition relates to dry material.

Properties

- Sugi absorbent materials are manu-factured using a high-grade material, based on cellulose.
- The swab material undergoes soft/elastic expansion immediately upon contacting aqueous solutions.
- Sugi absorbent materials absorb rapidly.
- Sugi absorbent materials will not adhere to wounds, which makes it easy to extract.
- The wound surface remains clean after the tampon has been removed.
- The absorptive capacity is approximately 17 times its own weight.

Range of Applications

- Postoperative nasal tampon and fixation splint for use after nasal surgery, especially after;
 - Plastic surgery of the septum
 - Septal rhinoplasty
 - Conchotomy
 - Paranasal sinus operations
- Tamponing after endoscopic ENT operations.
- Carrier medium for aqueous antibiotic or vasoconstrictive liquids.
- Fast acting tampon for use after severe nasal haemorrhaging.

Contraindications and Incompatibility

None known.

Application

The tampon can be cut to shape as required when dry. We recommend moistening the tampon with a small amount of physiological saline solution, aqueous antibiotic solution or vasoconstrictive solution prior to application.

Please note

- Sugi absorbent materials should be disposed of like bandaging materials.
- Sterility is only then guaranteed if the packaging is undamaged.
- Sugi absorbent materials should be stored in a dry place, at room temperature.
- Sugi Rhino absorbent swabs are radiolucent.
- Sugi absorbent materials are for use by qualified medical personnel only.
- In rare instances the physicochemical conditions associated with nasal surgery, both with and without nasal packing, may present a risk of toxic shock syndrome (TSS). Warning signs of TSS include: sudden fever, vomiting, diarrhea, dizziness, fainting (or near fainting when standing up), and/or a rash that looks like a sunburn. Should one of these symptoms occur, administer the appropriate medical treatment.

The information provided for Kettenbach products is based on comprehensive research and experience in application technology. Results are furnished to the best of our knowledge, subject to technical changes within the framework of product development. However, users must comply with and consider all recommendations and information in connection with any use.